

# WiseMini®

## Molecular Sieve Desiccant Strip




Two types:



Maximum length difference between two sachets in a strip roll is 1.5mm



Ultrasonic Seal is very strong and safe, which ensures each bag is cut at the right length when used in a dispensing machine

**GMP**

Our GMP Process flow ensures that all desiccant bags have equal moisture absorption properties on each reel



Full inspection product, promised for high quality



Molecular Sieve is able to absorb moisture at ultra-low humidity rates, which is great for most sensitive products



**WISESORBENT TECHNOLOGY LLC**

11 E Stow Rd, Marlton, NJ 08053, USA  
info@wisesorbent.com

WiseMini® molecular sieve desiccant strips eliminate the tedious process of manual insertion which is integral in high volume settings. Save on labor costs while simultaneously increase output via automatic dispensing. Meanwhile, our desiccants are completely safe to be in physical contact with any diagnostic products. Molecular Sieve is also essential for eliminating moisture at low humidity rates. Our desiccant strips are available in sizes from 0.25g to 10g while our dispensing and cutting machines can dispense up to 150 of our desiccant bags per minute.

Ultrasonic cut:



Hole cut:



### WISEMINI® MOLECULAR SIEVE DESICCANT STRIP STANDARD

ITEM		STANDARD
Loss on Drying ( 575±25°C )		≤ 4.5%
Moisture Absorption Capacity ( 77 °F /25°C )	RH = 40%	≥ 15.0%
	RH = 80%	≥ 16.5%

### WISEMINI® MOLECULAR SIEVE DESICCANT STRIP SPECIFICATION


Ultrasonic cut:

Spec	Product Code	W*L (mm\in)	Packaging Quantity
0.25g	1MSTY100A	16*30 \ 0.63*1.18	15000/ctn
0.5g	1MSTY101A	16*33.5 \ 0.63*1.32	15000/ctn
1g	1MSTY102A	19*42 \ 0.75*1.65	10000/ctn
2g	1MSTY103A	22*51.5 \ 0.87*2.03	5000/ctn
3g	1MSTY104A	26*60 \ 1.02*2.36	3500/ctn

Hole cut:

Spec	Product Code	W*L (mm\in)	Packaging Quantity
0.25g	2MSTG200A	17*30 \ 0.67*1.18	16000/ctn
0.5g	2MSTG201A	17*34 \ 0.67*1.34	15000/ctn
1g	2MSTG202A	22*40 \ 0.87*1.57	10000/ctn
2g	2MSTG203A	25*50 \ 0.98*1.97	5000/ctn
3g	2MSTG204A	25*60 \ 0.98*2.36	3500/ctn
5g	2MSTG205A	25*75 \ 0.98*2.95	2000/ctn
8g	2MSTG206A	37*75 \ 1.46*2.95	1200/ctn
10g	2MSTG207A	37*88 \ 1.46*3.46	1000/ctn

\* Container: 



• For additional desiccant solutions visit [wisesorbent.com](http://wisesorbent.com) •

**WISESORBENT TECHNOLOGY LLC**

11 E Stow Rd, Marlton, NJ 08053, USA  
info@wisesorbent.com

**Call today**  
to place an order  
**+1 (800) 272-5238**