


# WiseMop® Liquids Absorber



Strong water absorption (can absorb up to 350 times its own weight of distilled water, more than 40 times its own weight of normal saline)



Fast water absorption (can absorb 40-45 times its own weight in liquid within 35 seconds)  
\*The water absorption capacity is related to the concentration of metal salts



Water-soluble packaging, dissolved in water. Can be put into use directly without open the package



**WISORBENT TECHNOLOGY LLC**

11 E Stow Rd, Marlton, NJ 08053, USA  
info@wisesorbent.com

Used in hospitals, laboratories, nursing/elder care and other scenarios to quickly absorb/fix leaking liquids and prevent liquid contamination. It can absorb liquids within 1 min. After absorption, gel formation is easily to be collected and discard.



## WISEMOP® LIQUIDS ABSORBER STANDARD

ITEM	STANDARD
Loss on Drying	≤ 5.0%
Water absorption capacity	350-450 g/g
Absorption capacity for 0.9% NaCl solution	40-45 g/g
Absorption speed for the entire sachet	≤ 35s

## WISEMOP® LIQUIDS ABSORBER SPECIFICATION

\* Container: 

Spec	Product Code	W*L (mm\in)	Packaging Quantity
1g	2SPSR001A	25*42 \ 0.98*1.65	5000/5 gal tin
2g	2SPSR002A	25*52 \ 0.98*2.05	3000/5 gal tin
3g	2SPSR003A	35*55 \ 1.38*2.17	2500/5 gal tin
5g	2SPSR004A	35*65 \ 1.38*2.56	1500/5 gal tin



• For additional desiccant solutions visit [wisesorbent.com](http://wisesorbent.com) •

**WISORBENT TECHNOLOGY LLC**

11 E Stow Rd, Marlton, NJ 08053, USA  
info@wisesorbent.com

**Call today**  
to place an order  
**+1 (800) 272-5238**